

FDA Public Health Web Notification* : Complications Related to the Use of Bone Cement in Treating Compression Fractures of the Spine

Over the last decade, several types of complications have been reported in literature and to FDA following the use of polymethylmethacrylate (acrylic) bone cements to treat compression fractures of the spine. Labeling of marketed bone cements at the present time does not address use of this product to treat spinal compression fractures, nor does it address possible complications associated with this use. This document is intended to notify you of the types of reported complications associated with this use of bone cement.

Reported complications include serious injuries and deaths. Complications that relate specifically to the leakage of bone cements include soft tissue damage and nerve root pain and compression. Other reported complications generally associated with the use of bone cements in the spine include pulmonary embolism, respiratory and cardiac failure, abdominal intrusions/ileus, and death.¹⁻¹⁰ Each of these types of complications has been reported in conjunction with use of bone cements in both vertebroplasty and kyphoplasty procedures, two different techniques that employ bone cement to treat spinal compression fractures. Current information is insufficient to establish overall complication rates for either of these procedures or to distinguish the types and incidence of complications associated with the use of bone cements in each type of procedure.

Current Status

In the absence of labeling addressing use of bone cements to treat compression fractures of the spine, we encourage physicians to be aware of considerations and recommendations described in the literature and by professional organizations regarding patient selection, treatment techniques, potential complications, and patient monitoring when considering the use of bone cements in procedures to treat such fractures.^{5,11-15}

Because more study would be beneficial, we are currently working with professional organizations and manufacturers of orthopedic devices to develop a basis for evaluating the safety and effectiveness of bone cements used to treat compression fractures of the spine. Some manufacturers have now received approval to conduct clinical studies to address the safety and effectiveness of bone cements used in this way. Until such studies are complete and results are submitted to FDA, information available to us is too limited to perform, and we do not intend this notice to be, an assessment of the safety and effectiveness of bone cements used to treat compression fractures of the spine, the risks and benefits of this use of bone cement, or the surgeries that involve it.

Background

According to the National Osteoporosis Foundation, about 700,000 vertebral fractures occur annually; approximately 270,000 of these fractures are painful and clinically

* CDRH Web Notifications are intended to augment the existing Safety Notification program by providing a mechanism to quickly disseminate device-related safety information that may be beneficial to healthcare providers. Unlike other forms of Notifications, such as Safety Alerts or Public Health Advisories, Web Notifications usually do not make specific recommendations and are typically used in situations where the available information and our understanding of an issue are still evolving.

diagnosed. Most patients are treated non-operatively, but some do not respond to conservative treatment and are left with persistent pain and limited mobility. These patients are potential candidates for vertebroplasty or kyphoplasty procedures, two invasive procedures that use bone cement to treat spinal compression fractures.

There are differences in technique between vertebroplasty and kyphoplasty procedures, as well as among practitioners of each procedure, including differences in the selection, modification, and application of bone cements.

Vertebroplasty, developed in the 1980s, involves the percutaneous injection of a mixture of polymethylmethacrylate bone cement and a contrast agent, typically barium sulfate, into the vertebral bodies using fluoroscopic and/or computed tomography guidance.¹⁻¹⁷ Early vertebroplasty procedures were primarily performed to alleviate pain and to stabilize fractured bone in patients with hemangiomas, malignant metastases, or other types of tumors of the spine.

Kyphoplasty, developed in the 1990s, involves introducing a surgical instrument into the vertebral body, with the intent to elevate or expand the vertebra.^{3, 9, 18-20} Once this instrument is withdrawn, the space created is then filled with the bone cement mixture. By reducing and fixing the fracture in this way, kyphoplasty procedures may correct deformity and/or restore body height.

Acrylic bone cements have been used for many years for the fixation of metal and plastic prostheses in joint replacement and less frequently in the fixation of pathological fractures. These cements, however, have not been specifically evaluated for the treatment of spinal compression fractures. Existing cements designed for other uses are generally modified for use in treating spinal compression fractures. Modifications to cements may vary from physician to physician and among procedures. These types of modifications may include increasing the amount of contrast agent to improve x-ray visualization and changing the consistency and handling properties to address procedural goals (including changing the method of preparing the cement without altering its ingredients). To date, there are no standardized formulations, biomechanical standards or safety guidelines for the types or amounts of radio-opaque or other additives used, or for methods of preparing bone cement for use in the spine. The effects of modified bone cements on the spine and surrounding soft tissues have not been adequately studied to support marketing applications.

Reporting Adverse Events to FDA

The Safe Medical Devices Act of 1990 (SMDA) requires hospitals and other user facilities to report deaths and serious injuries associated with the use of medical devices, including bone cement. We request that you follow the procedures established by your facility for such mandatory reporting.

We also encourage you to report bone cement malfunctions. You can report these directly to the device manufacturer. You can also report to MedWatch, the FDA's

voluntary reporting program. You may submit reports to MedWatch one of four ways: online at <http://www.accessdata.fda.gov/scripts/medwatch/>; by telephone at 1-800-FDA-1088; by FAX at 1-800-FDA-0178; or by mail to MedWatch, Food and Drug Administration, HF-2, 5600 Fishers Lane, Rockville, MD 20857.

FDA Contact

Laura Alonge
Office of Surveillance and Biometrics (HFZ-510)
1350 Piccard Drive, Rockville, Maryland, 20850
Fax at 301-594-2968, or by e-mail at phann@cdrh.fda.gov.

Additionally, a voice mail message may be left at 301-594-0650 and your call will be returned as soon as possible.

References

1. Cortet B, Cotton A, Boutry N, et al. Percutaneous Vertebroplasty in the Treatment of Osteoporotic Vertebral Compression Fractures: An Open Prospective Study. *J Rheumatol*. 1999 Oct;26(1):2222-8.
2. Cyteval C, Sarrahere MP, Roux JO, et al. Acute osteoporotic vertebral collapse: open study on percutaneous injection of acrylic cement in 20 patients. *AJR* 1999;173(6):1685-90.
3. Garfin SR, Yuan HA, Reiley MA. New technologies in spine: kyphoplasty and vertebroplasty for the treatment of painful osteoporotic compression fractures. *Spine* 2001; 26(14):1511-5.
4. Grados F, Depriester C, Cayrolle G, et al. Long-term observations of vertebral osteoporotic fractures treated by percutaneous vertebroplasty. *Rheumatology* 2000; Dec;39(12):1410-14.
5. Heini PF, Walchli B, Berlemann U. Percutaneous transpedicular vertebroplasty with PMMA: operative technique and early results. A prospective study for the treatment of osteoporotic compression fractures. *Eur Spine J*. 2000 Oct;9(5):445-50.
6. Levine SA, Perin LA, HayesD, et al. An evidenced-based evaluation of percutaneous vertebroplasty. *Managed Care* 2000 Mar; 9(3):56-60, 63.
7. Martin JB, Jean B, Sugi K, et al. Vertebroplasty: clinical experience and follow-up results. *Bone* 1999 Aug;25(2 suppl):11S-15S.

8. Ratliff J, Nguyen T, Heiss J. Root and spinal cord compression from methylmethacrylate vertebroplasty. *Spine* 2001; 26(13):E300-2.
9. Watts NB, Harris ST, Genant HK. Treatment of painful osteoporotic vertebral fractures with percutaneous vertebroplasty or kyphoplasty. *Osteoporos Int.* 2001;12(6):429-37.
10. Padovani B, Kasriel O, Brunner P, et al. Pulmonary Embolism caused by acrylic cement: a rare complication of percutaneous vertebroplasty. *AJNR* 1999;20(3):375-7.
11. Do HM. Resonance imaging in the evaluation of patients for percutaneous vertebroplasty. *Top Magn Reson Imaging* 2000 Aug;11(4):235-44.
12. Hardouin P, Grados F, Cotton A, et al. Should percutaneous vertebroplasty be used to treat osteoporotic fractures? An update. *Joint Bone Spine* 2001 May;68(3):216-21.
13. Mathis JM, Barr JD, Belkoff SM, et al. Percutaneous vertebroplasty: a developing standard of care for vertebral compression fractures. *AJNR* 2001;22(2):373-81.
14. Contemporary Concepts in Spine Care – Percutaneous Vertebral Augmentation. North American Spine Society. Lagrange, Illinois. Dec. 2002.
15. Standard for the Performance of Percutaneous Vertebroplasty. American College of Radiology. Reston, Virginia. January 2001.
16. McGraw JK, Lippert JA, Minkus KD, et al. Prospective evaluation of pain relief in 100 patients undergoing percutaneous vertebroplasty results and follow-up. *JVIR* 2002;13:883-6.
17. Zoarski GH, Snow P, Olan WJ, et al. Percutaneous vertebroplasty for osteoporotic compression fractures: quantitative prospective evaluation of long-term outcomes. *J Vasc Interv Radiol* 2002;13:139-48.
18. Hardouin P, Fayada P, Leclet H, et al. Kyphoplasty. *Joint Bone Spine* 2002 May;69(3):256-61.
19. Lieberman IH, Dudeney S, Reinhardt MK, et al. Initial outcome and efficacy of kyphoplasty in the treatment of painful osteoporotic vertebral compression fractures. *Spine* 2001;26(14):1631-7.
20. Theodorou DJ, Theodorou SJ, Duncan TD, et al. Percutaneous balloon kyphoplasty for the correction of spinal deformity in painful vertebral compression fractures. *Clin Imaging* 2002 Jan-Feb;26(1):1-5.